


Liquid Bath Accessories



Mechanical support for probes

When setting up a new calibration bath, you need a way to suspend your probes in the bath fluid. We recommend our modular mechanical support systems. Made of fine-quality steel and machined parts, these support components combine in hundreds of ways to solve almost any probe suspension problem.

Single probe kit

Our single probe kit is a good way to get started. It has one medium clamp, one 10-inch rod, one bosshead, and one V-base for holding one probe.

Economy kit

Our economy kit includes two V-bases, one 29-inch stainless steel rod, two 23-inch rods, five bossheads, two micro-clamps, and one medium clamp. This will build a bench-mounted frame for suspending one rod and three clamps over a bath opening. Simply add to your setup as needed. Choose from any of the listed accessories.



Individual hardware

Our selection of clamps and stands provides a simple way to hold probes and thermometers in baths during calibration. System components can be assembled in a number of ways to suit individual needs. Select clamps, rods, bossheads, and bases to fit your individual needs.

Ordering Information		
Selected kits		Model
Single probe kit	Includes: 1 medium clamp, 1 10-inch rod, 1 bosshead, and 1 V-base	2051
Economy kit	Includes: 2 V-bases, 1 29-inch stainless steel rod, 2 23-inch rods, 5 bossheads, 2 micro-clamps, and 1 medium clamp	2050
Individual hardware		
Micro-clamps	Holds thermometers and probes with diameters to 0.75 inch.	2055 (Pkg. of 2)
Medium clamp	Holds diameters up to 1.75 inches.	2056 (Pkg. of 1)
Non-slip tape	Increases grip on clamps and bossheads.	2057 (1 roll)
Stainless steel rods	Used to assemble supports, frameworks, or scaffolds.	2058 10" (Pkg. of 1)
		2059 20" (Pkg. of 1)
		2060 23" (Pkg. of 1)
		2061 29" (Pkg. of 1)
Bossheads	Clamps two rods at right angles. Also attaches clamps to rods.	2062 (Pkg. of 5)
Screw bases	Holds one rod and is screwed to the surface of your bench or bath lid. 2.5-inch-diameter base; screws included.	2063 (Pkg. of 4)
V-base	Holds one rod. Weighted for excellent stability. 2.2 pounds (1 kg).	2064 (Pkg. of 1)
Large V-base	Same as above but larger. Recommended for holding SPRTs and large probes. 4.4 pounds (2 kg). (May be too large for some baths.)	

